

产品分析证书
Certificate of Analysis

产品信息(Product Information)

产品编号 (CAT No.): N0107

中文名称: N-反式-对-香豆酰基去甲辛弗林


English Name: Octopamine; N-p-coumaroyl-

CAS No.: 66648-45-1

分子式 (M. F.): C₁₇H₁₇NO₄

分子量 (M. W.): 299.326

批号 (Batch No.): N0107-1g-M230522101833



检测项目及结果 (Specifications&results)

检验项目 (Test Item)	检验规格 (Specifications)	检验结果 (Results)
颜色 Color	白色至类白色 (White to Off-white)	类白色 (Off-white)
外观 Appearance	粉末或固体 (Powder or Solid)	粉末 (Powder)
纯度 Purity	98%	Conforms
质谱 Mass	Comply with the structure	Conforms
核磁 NMR	Comply with the structure	Conforms

贮存条件 (Storage): 2~8℃, protected from light, keep package airproofed. <-15℃ when not in use for a long time. (密封保存, 置于阴凉干燥处, 长期不使用最好放置于-15℃以下环境中)

温馨提示 (Reminder): The sample solutions should be prepared immediately before use. (请现配现用)



QA: Joy Li

产品售后 (After-sale service): In case of quality issue, please contact us within 15 days after receipt of the product. (若对于我们的产品质量有任何疑问, 请于收到货后 15 天内联系我们)

(We appreciate you much when publishing your paper and marking the chemical purchase from Naturewill Biotechnology Co., Ltd. (Chengdu, China) 当发表论文时, 请您标注化合物购买于四川恒诚致远生物科技有限公司, 非常感谢您的支持。)

样品名称:104276-230501N-对香豆酰真蛸胺

仪器方法:N-对香豆酰真蛸胺

稀释倍数:1

报告方法:---

进样时间:2023/5/18 13:35:28

数据处理时间:2023/5/18 15:21:05

进样功能:进样

样品瓶:Left.A4

进样体积:5μL

样品类型:未知样


处理方法:默认处理方法

处理者:admin

进样结束时间:2023/5/18 14:06:11

序列名称:N-对香豆酰真蛸胺-2023/5/18 10:31:42

谱图



峰结果表					
序号	保留时间	峰面积	峰面积%	峰高	理论塔板数
1	10.822	2.931	0.034	0.246	22077
2	13.939	8590.491	99.694	563.036	19463
3	17.015	2.209	0.026	0.132	20232
4	22.829	19.278	0.224	0.778	18753
5	25.949	1.989	0.023	0.100	23262
统计	90.554	8616.899	100.000	564.291	

