


SAMPLE INFORMATION

Sample Name:	Ursolic acid	Acquired By:	hc01
Sample Type:	Unknown	Sample Set Name:	
Vial:	15	Acq. Method Set:	Ursolic acid
Injection #:	1	Processing Method	sample
Injection Volume:	10.00 ul	Channel Name:	SATIN
Run Time:	25.0 Minutes	Proc. Chnl. Descr.:	
Date Acquired:	2025-5-6 10:22:28 CST		
Date Processed:	2025-5-6 10:49:03 CST		

Auto-Scaled Chromatogram



Peak Results

	Name	RT	Area	% Area	USP Plate Count	USP Resolution
1		9.269	10718985	99.94	13510.63	
2		12.022	1947	0.02	561425.07	14.17
3		12.885	2980	0.03	615908.91	12.68
4		19.934	1677	0.02	1488741.02	102.22