

产品分析证书 Certificate of Analysis

产品信息(Product Information)

产品编号 (CATNo.): N0852

中文名称: 人参炔醇

English Name: Panaxynol

CAS No.: 81203-57-8

分子式 (M. F.): C₁₇H₂₄O

分子量 (M. W.): 244.378

批号 (Batch No.): N0852-10mg-M231213112557



检测项目及结果 (Specifications&results)

检验项目 (Test Item)	检验规格 (Specifications)	检验结果 (Results)
颜色 Color	淡黄色至黄色 (light yellow to Yellow)	淡黄色 (light yellow)
外观 Appearance	油状液体 (oil liquid)	油状液体 (oil liquid)
纯度 Purity	98%	Conforms
质谱 Mass	Comply with the structure	Conforms
核磁 NMR	Comply with the structure	Conforms

贮存条件 (Storage): 2~8°C, protected from light, keep package airtight. <-15°C when not in use for a long time. (密封保存, 置于阴凉干燥处, 长期不使用最好放置于-15°C以下环境中)


温馨提示 (Reminder): The sample solutions should be prepared immediately before use. (请现配现用)



QA: Chao Zhang

产品售后 (After-sale service): In case of quality issue, please contact us within 15 days after receipt of the product. (若对于我们的产品质量有任何疑问, 请于收到货后 15 天内联系我们)

(We appreciate you much when publishing your paper and marking the chemical purchase from Naturewill Biotechnology Co., Ltd. (Chengdu, China) 当发表论文时, 请您标注化合物购买于四川恒诚致远生物科技有限公司, 非常感谢您的支持。)



峰结果

名称	保留时间 (分钟)	面积 (微伏*秒)	高度 (微伏)	% 面积
1 峰1	5.992	193001	22505	1.28
2 峰2	8.173	572249	53618	3.80
3 峰3	16.377	123278	7110	0.82
4 峰4	19.499	398737	19212	2.65
5 峰5	25.209	209547	13182	1.39
6 峰6	25.826	329207	20368	2.19
7 峰7	32.655	87778	6342	0.58
8 峰8	33.301	105441	6935	0.70
9 峰9	33.757	618591	48913	4.11
10 峰10	34.090	11099280	580277	73.70
11 峰11	35.108	1084036	64287	7.20
12 峰12	35.835	134628	8539	0.89