

产品分析证书

Certificate of Analysis

产品信息(Product Information)

产品编号 (CAT No.) : N1822

中文名称: 刺五加苷 E


English Name: Eleutheroside E

CAS No.: 39432-56-9

分子式 (M. F.): C₃₄H₄₆O₁₈

分子量 (M. W.): 742.72

批号 (Batch No.): 2207182201



检测项目及结果 (Specifications&results)

检验项目 (Test Item)	检验规格 (Specifications)	检验结果 (Results)
颜色 Color	类白色至淡黄色(Off-white to light yellow)	类白色(Off-white)
外观 Appearance	粉末或固体(Powder or Solid)	粉末(Powder)
纯度 Purity	≥98.0% by HPLC	Conforms
质谱 Mass	Comply with the structure	Conforms
核磁 NMR	Comply with the structure	Conforms

贮存条件 (Storage) : 2~8℃, protected from light, keep package airproofed. <-15℃ when not in use for a long time. (密封保存, 置于阴凉干燥处, 长期不使用最好放置于-15℃以下环境中)

温馨提示 (Reminder) : The sample solutions should be prepared immediately before use.(请现配现用)


QC: Ben Yang

QA: Joy Li



产品售后 (After-sale service): In case of quality issue, please contact us within 15 days after receipt of the product. (若对于我们的产品质量有任何疑问, 请于收到货后 15 天内联系我们)

HPLC Chromatogram



峰序	组分名	保留时间[min]	峰高[uV]	峰面积[uV*s]	峰面积% 含量[%]	峰类型
1		2.457	41	143	0.01000	0.00000 BB
2		4.337	499	5997	0.42068	0.00000 VP
3		9.387	35	95	0.00663	0.00000 VV
4		16.738	24	106	0.00743	0.00000 BB
5		24.813	45501	1418563	99.51092	0.00000 VV
6		26.327	62	632	0.04434	0.00000 VV