

产品分析证书 Certificate of Analysis

产品信息(Product Information)

产品编号 (CAT No.): N0012

中文名称: 单咖啡酰酒石酸


English Name: Caftaric acid

CAS No.: 67879-58-7

分子式 (M. F.): C₁₃H₁₂O₉

分子量 (M. W.): 312.23

批号 (Batch No.): N2500120321



检测项目及结果 (Specifications&results)

检验项目(Test Item)	检验规格 (Specifications)	检验结果 (Results)
颜色 Color	淡黄色至黄色(light yellow to Yellow)	淡黄色(light yellow)
外观 Appearance	粉末或固体(Powder or Solid)	粉末(Powder)
纯度 Purity	98%	Conforms
质谱 Mass	Comply with the structure	Conforms
核磁 NMR	Comply with the structure	Conforms

检测信息 (Detection Information)

色谱柱(Column)	Hypersil ODS-2 250*4.6mm 5μm
流动相(Mobile phase)	甲醇: 0.1% 磷酸= 13-20% 7min-30% 18min 甲醇 (Methanol : 0.1% Phosphoric acid in water =13-20% 7min-30% 18min Methanol)
检测波长(Detection mode)	327nm
柱温 (Column temperature)	30℃

贮存条件 (Storage): 2~8℃, protected from light, keep package airproofed. <-15℃ when not in use for a long time. (密封保存, 置于阴凉干燥处, 长期不使用最好放置于-15℃以下环境中)

温馨提示 (Reminder): The sample solutions should be prepared immediately before use.(请现配现用)




QA: Joy Li

产品售后 (After-sale service): In case of quality issue, please contact us within 15 days after receipt of the product. (若对于我们的产品质量有任何疑问, 请于收到货后 15 天内联系我们)

(We appreciate you much when publishing your paper and marking the chemical purchase from Naturewill Biotechnology Co., Ltd.(Chengdu,China) 当发表论文时, 请您标注化合物购买于四川恒诚致远生物科技有限公司, 非常感谢您的支持。)

HPLC



处理通道说明: MWD: 信号 A, 327 nm/带宽: 4 nm

	处理通道说明	保留时间 (分钟)	面积 (微伏·秒)	% 面积	高度 (微伏)
1	MWD: 信号 A, 327 nm/带宽 4 nm	10.315	3398	0.02	360
2	MWD: 信号 A, 327 nm/带宽 4 nm	12.266	15066373	99.37	1421495
3	MWD: 信号 A, 327 nm/带宽 4 nm	13.483	2193	0.01	257
4	MWD: 信号 A, 327 nm/带宽 4 nm	16.480	11493	0.08	1121
5	MWD: 信号 A, 327 nm/带宽 4 nm	16.939	3665	0.02	401
6	MWD: 信号 A, 327 nm/带宽 4 nm	18.550	5567	0.04	520
7	MWD: 信号 A, 327 nm/带宽 4 nm	19.534	7956	0.05	367
8	MWD: 信号 A, 327 nm/带宽 4 nm	23.425	11093	0.07	1016
9	MWD: 信号 A, 327 nm/带宽 4 nm	25.278	50189	0.33	4626

NMR

