

产品分析证书 Certificate of Analysis

产品信息(Product Information)

产品编号 (CAT No.): N0110

中文名称: 柚皮苷


English Name: Naringin

CAS No.: 10236-47-2

分子式 (M. F.): C₂₇H₃₂O₁₄

分子量 (M. W.): 580.539

批号 (Batch No.): N2401011022



检测项目及结果 (Specifications&results)

检验项目 (Test Item)	检验规格 (Specifications)	检验结果 (Results)
颜色 Color	白色至类白色 (White to Off-white)	白色 (White)
外观 Appearance	粉末或固体 (Powder or Solid)	固体 (Solid)
纯度 Purity	98%	Conforms
质谱 Mass	Comply with the structure	Conforms
核磁 NMR	Comply with the structure	Conforms

贮存条件 (Storage): 2~8°C, protected from light, keep package airtight. <-15°C when not in use for a long time. (密封保存, 置于阴凉干燥处, 长期不使用最好放置于-15°C以下环境中)


温馨提示 (Reminder): The sample solutions should be prepared immediately before use. (请现配现用)



QA: Chao Zhang

产品售后 (After-sale service): In case of quality issue, please contact us within 15 days after receipt of the product. (若对于我们的产品质量有任何疑问, 请于收到货后 15 天内联系我们)

(We appreciate you much when publishing your paper and marking the chemical purchase from Naturewill Biotechnology Co., Ltd. (Chengdu, China) 当发表论文时, 请您标注化合物购买于四川恒诚致远生物科技有限公司, 非常感谢您的支持。)



=====
 Area Percent Report
 =====

Sorted By : Signal
 Multiplier : 1.0000
 Dilution : 1.0000
 Use Multiplier & Dilution Factor with ISTDs


Signal 1: VWD1 A, Wavelength=220 nm

Peak #	RetTime [min]	Type	Width [min]	Area [mAU*s]	Height [mAU]	Area %
1	7.469	MM	0.1107	8.95209	1.34792	0.1050
2	7.743	MM	0.1071	25.79180	4.01219	0.3026
3	8.230	MM	0.0765	8440.45703	1839.63306	99.0416
4	8.560	MF	0.0706	32.78051	7.73841	0.3847

MassPeaks:467

Spectrum Mode:Averaged 0.945-0.998(75-79) Base Peak:581.4(410899)

BG Mode:Calc Segment 1 - Event 1



MassPeaks:435

Spectrum Mode:Averaged 0.958-1.012(76-80) Base Peak:579.4(1059555)

BG Mode:Calc Segment 1 - Event 2

