

产品分析证书
Certificate of Analysis

产品信息(Product Information)

产品编号 (CAT No.): N3168


中文名称: 桔霉素

CAS No.: 518-75-2

分子式 (M. F.): C₁₃H₁₄O₅

分子量 (M. W.): 250.25

批号 (Batch No.): 2403316801



检测项目及结果 (Specifications&results)

检验项目 (Test Item)	检验规格 (Specifications)	检验结果 (Results)
颜色 Color	Yellow crystalline powder	Complies
外观 Appearance	粉末或固体 (Powder or Solid)	粉末 (Powder)
纯度 Purity	98%	Conforms
质谱 Mass	Comply with the structure	Conforms
核磁 NMR	Comply with the structure	Conforms

贮存条件 (Storage): 2~8℃, protected from light, keep package airproofed. <-15℃ when not in use for a long time. (密封保存, 置于阴凉干燥处, 长期不使用最好放置于-15℃以下环境中)

温馨提示 (Reminder): The sample solutions should be prepared immediately before use. (请现配现用)




QA: Joy Li

产品售后 (After-sale service): In case of quality issue, please contact us within 15 days after receipt of the product. (若对于我们的产品质量有任何疑问, 请于收到货后 15 天内联系我们)

(We appreciate you much when publishing your paper and marking the chemical purchase from Naturewill Biotechnology Co., Ltd. (Chengdu, China) 当发表论文时, 请您标注化合物购买于四川恒诚致远生物科技有限公司, 非常感谢您的支持。)

HPLC

Additional Info : Peak(s) manually integrated



Area Percent Report

Sorted By : Signal
Multiplier: : 1.0000
Dilution: : 1.0000
Use Multiplier & Dilution Factor with ISTDs

Signal 1: DAD1 B, Sig=340,4 Ref=360,100

Peak #	RetTime [min]	Type	Width [min]	Area [mAU*s]	Height [mAU]	Area %
1	6.722	BB	0.1060	4.72117	7.05950e-1	0.3569
2	8.555	BB	0.1326	1317.96423	152.80836	99.6431

Totals : 1322.68541 153.51431

2023020621 299 (11.044)

1: Scan ES+
1.89e6

