


SAMPLE INFORMATION

Sample Name: Wogonin	Acquired By: System
Sample Type: Unknown	Sample Set Name:
Vial: 54	Acq. Method Set: Wogonin
Injection #: 1	Processing Method: Samples
Injection Volume: 10.00 ul	Channel Name: 274.0nm
Run Time: 28.0 Minutes	Proc. Chnl. Descr.: PDA 274.0 nm
Date Acquired: 2022-8-22 12:23:33 CST	
Date Processed: 2022-8-22 13:38:46 CST	

Auto-Scaled Chromatogram



Peak Results

	RT	Area	% Area	USP Plate Count	USP Resolution
1	1.786	62756	0.16		
2	2.358	60915	0.16		
3	6.532	22290	0.06		
4	7.339	146107	0.38		
5	7.835	17611	0.05		
6	11.046	59396	0.15		
7	11.572	38475188	99.05		