


产品信息(Product Information)

批号 (Batch No.): N2407087633



检测项目及结果 (Specifications&results)

检验项目(Test Item)	检验规格 (Specifications)	检验结果 (Results)
颜色 Color	Purple brown	Purple brown
外观 Appearance	粉末或固体(Powder or Solid)	粉末(Powder)
纯度 Purity	98%	Conforms
质谱 Mass	Comply with the structure	Conforms
核磁 NMR	Comply with the structure	Conforms

检测信息(Detection Information)

色谱柱(Column)	Hypersil ODS-2 250*4.6mm 5μm
流动相(Mobile phase)	甲醇：0.1%磷酸水溶液=72:28
检测波长(Detection mode)	203nm
柱温 (Column temperature)	30℃

温馨提示 (Reminder): The sample solutions should be prepared immediately before use. (请现配现用)



QA: Chao Zhang

(We appreciate you much when publishing your paper and marking the chemical purchase from Naturewill Biotechnology Co., Ltd.(Chengdu,China) 当发表论文时, 请您标注化合物购买于四川恒诚致远生物科技有限公司, 非常感谢您的支持。)