

产品分析证书

Certificate of Analysis

产品信息(Product Information)

产品编号 (CAT No.): N0239

中文名称: 淫羊藿苷


English Name: Icariin

CAS No.: 489-32-7

分子式 (M. F.): C₃₃H₄₀O₁₅

分子量 (M. W.): 676.668

批号 (Batch No.): 2206023901



检测项目及结果 (Specifications&results)

检验项目 (Test Item)	检验规格 (Specifications)	检验结果 (Results)
颜色 Color	淡黄色至黄色 (light yellow to Yellow)	黄色 (yellow)
外观 Appearance	粉末或固体 (Powder or Solid)	粉末 (Powder)
纯度 Purity	98%	Conforms
质谱 Mass	Comply with the structure	Conforms
核磁 NMR	Comply with the structure	Conforms

检测信息 (Detection Information)

色谱柱 (Column)	Aglient 250mm
流动相 (Mobile phase)	0.1%磷酸、乙腈
检测波长 (Detection mode)	270nm
柱温 (Column temperature)	30℃

贮存条件 (Storage): 2~8℃, protected from light, keep package airproofed. <-15℃ when not in use for a long time. (密封保存, 置于阴凉干燥处, 长期不使用最好放置于-15℃以下环境中)

温馨提示 (Reminder): The sample solutions should be prepared immediately before use. (请现配现用)



QA: Joy Li


产品售后 (After-sale service): In case of quality issue, please contact us within 15 days after receipt of the product. (若对于我们的产品质量有任何疑问, 请于收到货后 15 天内联系我们)

(We appreciate you much when publishing your paper and marking the chemical purchase from Naturewill Biotechnology Co., Ltd. (Chengdu, China) 当发表论文时, 请您标注化合物购买于四川恒诚致远生物科技有限公司, 非常感谢您的支持。)

SAMPLE INFORMATION

Sample Name:	Icariin	Acquired By:	System
Sample Type:	Unknown	Sample Set Name:	
Vial:	62	Acq. Method Set:	Icariin
Injection #:	1	Processing Method:	Samples
Injection Volume:	10.00 ul	Channel Name:	270.0nm
Run Time:	25.0 Minutes	Proc. Chnl. Descr.:	PDA 270.0 nm
Date Acquired:	2023-11-18 10:27:29 CST		
Date Processed:	2023-11-18 11:10:08 CST		

Auto-Scaled Chromatogram



Peak Results

	Name	RT	Area	% Area	USP Plate Count	USP Resolution
1		8.374	15540	0.11	37753.11	
2		9.122	8170	0.06	38405.96	4.34
3		11.633	6154	0.04	60170.18	13.50
4		12.278	14364651	99.79	69862.18	3.33